

Application Note


Abstract

With the growing price of oil, many alternative energy sources are being explored. Natural gas prices are actually decreasing, in large part to the use of hydraulic fracturing in areas like the Marcellus Shale in Pennsylvania. Due to increased concern over the hydrofracturing process and the release of methane and other chemicals into the local drinking water, a need has developed for fast and accurate analysis of methane in water. This study will evaluate a method developed for the determination of methane, ethane, ethene, and propane in water using a Purge and Trap concentrator, autosampler, and GC/FID. Comparisons to current headspace methods, RSK 175¹ and BOL6019², will be made.




Introduction

Hydraulic fracturing, also known as “fracking”, is a drilling process currently used to recover natural gas from sources like coalbeds and shale formations. This involves injecting large amounts of water, mixed with sand and chemicals, at high pressures to break up the shale to release the gas. An outline of the hydraulic fracturing process can be found in Figure 1. Natural gas from these hard-to-extract sources is becoming increasingly popular and is projected to grow to nearly 45% of the nation’s natural gas supplies by 2035³.



Hydraulic Fracturing

A new way of drilling for natural gas



Courtesy of Ohio Department of Natural Resources Division of Mineral Resources Management (Illustration modified by TOMWC)


Figure 1: Process of Hydraulic Fracturing to Drill for Natural Gas⁴

Even though hydraulic fracturing is a relatively old practice, first employed over 60 years ago to drill for oil in Oklahoma⁵, there has been little research into the impact of its increasing use as a drilling process for natural gas. With growing concern over the environmental effects of fracking on water quality, the United States EPA has begun studies to monitor the treatment methods and environmental impact with the goal of standardization by 2013/2014⁶. Figure 2 shows a map of shale gas formations in the continental United States indicating the potential widespread environmental impact the fracking process could have.

The current method for determining natural gas constituents (methane, ethane, and ethene) in water is RSK 175¹. Propane has been added to this list in modified methods such as BOL6019², developed by the Pennsylvania Department of Environmental Protection (PADEP). These analyses require headspace gas chromatography. This study will also utilize a flame ionization detector (FID), although RSK 175¹ also allows thermal conductivity (TCD) as well as electron capture detectors (ECD) to be used. Due to the lack of EPA methods requiring headspace analysis, these instruments are not normally found in

Sales/Support: 800-874-2004 · **Main:** 513-229-7000
4736 Socialville Foster Rd., Mason, OH 45040
www.teledynetekmar.com

environmental laboratories. This poster demonstrates an alternative analysis using purge and trap concentration, which are typically available in many environmental labs.



Source: Energy Information Administration based on data from various published studies.
Updated: March 10, 2010

Figure 2: Shale Gas Formations in the United States⁷

For this study, a Stratum Purge and Trap Concentrator (PTC) was used in conjunction with an AQUATek 100 Autosampler. This set-up allows for complete automation of sample preparation for the analysis of liquid samples for purge and trap. A recirculating chiller bath was also required to maintain a sample temperature of less than 4°C. This technique also requires a 25mL purge volume.

Utilizing a GC/FID, a linear calibration was performed and percent Relative Standard Deviation (%RSD) and Method Detection Limits (MDLs) were determined for the full list of compounds. This purge and trap technique, which requires no sample manipulation unlike RSK 175¹ and BOL6019², will be evaluated using real world replicate samples versus the current headspace method. Similarly to BOL6019², calibrations will be performed on aqueous standards rather than the gaseous standards used in RSK 175¹.

Experimental-Instrument Conditions

The Stratum PTC and AQUATek 100 Autosampler were coupled to a GC/FID for analysis. A new proprietary trap was also utilized. The GC was configured with a Restek Rt-U-PLOT 30m x 0.53mm column. The GC/FID parameters are outlined in Tables 1 and 2. Table 3 outlines the P&T and autosampler conditions. Due to the extreme volatility of these analytes, a recirculating bath was employed to chill the samples below 4°C.

Sales/Support: 800-874-2004 · **Main:** 513-229-7000
4736 Socialville Foster Rd., Mason, OH 45040
www.teledynetekmar.com

GC Parameters	
GC:	GC/FID
Column:	Restek Rt-U-PLOT 30m x 0.53mm
Oven Program:	Initial Temp 50°C, 10°C/min to 190°C
Inlet:	250°C
Gas:	Helium
Pressure:	19.44psi
Split Ratio:	20:1

FID Parameters	
Temperature:	250°C
Hydrogen Flow:	40mL/min
Air Flow:	400mL/min
Mode:	Constant Makeup Flow
Makeup Flow:	30mL/min
Makeup Gas:	Helium

Tables 1 & 2: GC and FID Parameters

Stratum PTC and AQUATek 100 Parameters			
Variable	Value	Variable	Value
Pressurize Time	0.85 min	Purge Time	1.5 min
Sample Transfer Time	1.00 min	Purge Temp	20°C
Rinse Loop Time	0.85 min	Purge Flow	20mL/min
Sweep Needle Time	0.35 min	Dry Purge Time	0.0 min
Bake Rinse	On	Dry Purge Temp	20°C
Bake Rinse Cycles	1	Dry Purge Flow	100mL/min
Bake Rinse Drain Time	0.60 min	GC Start	Start of Desorb
Presweep Time	0.35 min	Desorb Preheat Temp	95°C
Water Temp	90°C	Desorb Drain	On
Valve Oven Temp	80°C	Desorb Time	2.00 min
Transfer Line Temp	80°C	Desorb Temp	100°C
Sample Mount Temp	60°C	Desorb Flow	300mL/min
Purge ready Temp	35°C	Bake Time	8.00 min
Condenser Ready Temp	40°C	Bake Temp	100°C
Condenser Purge Temp	20°C	Bake Flow	400mL/min
Standby Flow	10mL/min	Condenser Bake Temp	200°C
Pre-Purge Time	0.5 min		
Pre-Purge Flow	40.0mL/min		
Sample Heater	Off		
Sample Preheat Time	1.00 min		
Sample Temp	40°C		

Table 3: Stratum PTC and AQUATek 100 Parameters (Stratum PTC Parameters are in Blue)

Calibration Data

To make the stock solutions, a 500mL volumetric flask filled with de-ionized water was placed in an ice water bath and purged with a reference gas corresponding to each of the four analytes. Each gas was bubbled through chilled water for two hours to make individual concentrated standards. Unlike in RSK 175¹, calibrations in BOL6019² are performed using aqueous rather than gaseous standards. This study also employs an aqueous calibration, where standards are analyzed under the same conditions as samples.

Calibration standards were made from serial dilutions of these stock standards by calculating the concentration of saturated gas solutions in water at 0°C. For example, the saturation point of methane in 0°C water is 39.59 mg/L at atmospheric pressure. Calibration standards were made in 50mL volumetric

flasks filled to volume with chilled de-ionized water over a range of 8 ppb to 40 ppm. Samples were transferred to headspace free 40mL vials for analysis. These standards were additionally made for the remaining compounds of interest. Samples for the comparative headspace analysis were prepared similarly and added to 22mL headspace vials for analysis.

The calibration data generated during this study was evaluated by percent Relative Deviation (%RSD). Method detection limits were also established for all compounds by analyzing replicates lab fortified blanks (LFBs). Calibration data can be found in Table 4. In addition, an example of a chromatogram for a 40ppm methane standard can be found in Figure 3. A blank analyzed after the highest calibration standard was used to calculate the percent carryover which was less than 0.5% for all compounds.

Compound	Calibration Range	Relative Response Factor (RRF)	% Relative Deviation (%RSD)	Method Detection Limit (MDL)
Methane	20ppb to 40ppm	573	10.2	8.3 ppb
Ethane	26ppb to 66ppm	4024	11.1	10.5 ppb
Ethene	67ppb to 281ppm	2125	6.6	N/A
Propane	17ppb to 71ppm	8424	12.1	2.7 ppb

Table 4: Calibration and MDL Data for Methane, Ethane, Ethene, and Propane




Figure 3: Chromatogram of 40ppm Methane Standard

Real world drinking water samples were run head-to-head using both the current headspace method and the proposed purge and trap. A set of 37 replicate samples were analyzed using similar conditions for both methods. The average bias for these replicate samples for the purge and trap method was +4% for methane and -5.6% for ethane. Ethene and propane are rarely found in drinking water samples. These similar results show the quicker and more efficient purge and trap method to be an acceptable alternative to the current headspace method used at PADEP, BOL6019².

Conclusions

With increased interest in alternative energy sources, hydraulic fracturing has become a common practice in the extraction of natural gas from coalbeds and shale formations across the United States. Unfortunately, there has not been adequate time to measure the environmental impact of these procedures. Regulatory agencies are looking for easy and reliable testing methods to monitor these effects.

This study demonstrates a new method for analyzing these gases using the Teledyne Tekmar Stratum PTC and AQUATEk 100 Autosampler coupled with a GC/FID system. This method met all performance criteria outlined in the current headspace methods, RSK 175¹ and BOL6019². While this purge and trap method has been validated in this study, further method development and optimization are underway at Teledyne Tekmar. Additional comparisons to the existing headspace methods are also in progress, including calibration requirements, sample recovery, and supportive quality control data.

By completely automating the sample preparation in the purge and trap method, efficiency and throughput can be greatly increased while saving time and money. The cycle time for the purge and trap method is 15 minutes (four samples per hour), while the 30 minute equilibration time used in the headspace method pushes the cycle time out roughly 45 minutes (one sample per hour). There is no need to manipulate the samples required by the headspace method which eliminates the potential for human error and employs instrumentation many environmental laboratories already possess.

Teledyne Tekmar would like to thank the Pennsylvania Department of Environmental Protection for their participation in this study.

References

1. "Analysis of Dissolved Methane, Ethane, and Ethylene in Ground Water by a Standard Gas Chromatographic Technique", Don H. Kampbell and Steve A. Vandegrift, *Journal of Chromatographic Science*, Vol. 36, 253256, May 1998.
2. "C1 to C3 Dissolved Gases via Headspace and GC/FID, Revision 5", PA Dept. of Environmental Protection, Feb 2011.
3. "Annual Energy Outlook 2009", Energy Information Administration, [http://www.eia.gov/oiaf/archive/aeo09/pdf/0383\(2009\).pdf](http://www.eia.gov/oiaf/archive/aeo09/pdf/0383(2009).pdf)
4. "Natural Gas Drilling and Water", Tip of the Mitt Watershed Council, <http://www.watershedcouncil.org/learn/hydraulic-fracturing/>
5. "The Facts about Hydraulic Fracturing" Ohio Department of Natural Resources, <http://www.ohiodnr.com/Portals/11/pdf/fracking-fact-sheet.pdf>
6. Hydraulic Fracturing, United States Environmental Protection Agency, <http://www.epa.gov/hydraulicfracturing>
7. "Review of Emerging Resources: U.S. Shale Gas and Shale Oil Plays", Energy Information Administration, <http://www.eia.gov/analysis/studies/usshalegas/>