

Part Number	Description
<b>Spare Parts - Glassware &amp; Chemistry Area</b>	

14-9751-067	Loop, MD Dryer, Teflon
14-9536-024	UV Lamp, Fusion
14-7680-124	Sparge Tube
14-7033-043	O-ring, sparge tube, in UV reaction vessel
100-446	O-Ring, Lamp, UV Reaction Vessel
14-8809-000	Tool for UV Lamp Lead Pins
14-9791-079	Connector for UV Lamp Guard
14-7021-100	UV Reactor Cap Assembly (PEEK)
14-7022-080	Metal plate for UV reaction vessel cap
110-206	UV Reactor Screws (peek cap)
050-114	Plug Con 6 CKT Amp 1480704-0
	UV Lamp Accessory Kit - Contains 6 O-Rings for Lamp, 3 O-rings for the Sparge Tube, 6 wire ties, 4 reactor cap screws, and a Pin Extractor for the UV leads.
14-9961-000	Reaction Vessel Glassware for UV Reactor.
14-9556-024	Fusion Front Cover (Tinted)
14-9564-094	25mL Mist Trap
14-9558-024	U-Tube Scrubber Assembly
14-9758-000	U-Shaped Scrubber Tube for Corrosive Scrubber.
14-9756-079	Adapter, 5/16-24 Male to 1/4-28 Female
14-9757-016	O-ring, 3-stage needle
14-8254-043	Copper-20 mesh, 1.5 oz Kit
511-895	Tin 20 mesh 1.5 oz
511-876	Pyrex wool 1 gram
511-738	Sparger, 25mL, 1" Neck
14-9884-024	O-ring .734ID x .103 Wall Fluoro Viton
14-9893-043	Seal, O-ring, .206" ID x .103" W Viton
14-2007-043	Fitting, 1" Tube - 23 mL Sparger
14-9885-016	Fitting, 1/2" Tube - 25 mL Sparger, Kynar
14-9886-016	SO3 Scrubber Kit, Fusion
14-8203-200	Kit, Mist Catcher

<b>Spare Parts - Pumper &amp; Access.</b>	
---	--

14-9719-000	Pumper, Kloehn, V6
14-9784-050	7-port valve for syringe pumper
14-7127-152	Syringe, 25ml with plunger
14-8661-001	Screw, 6-32 x 1-5 SHCS
14-7127-000	Pin, Drive, for 14-7127-052 Syringe

<b>Spare Parts - Detector and Electronics</b>	
---	--

14-9561-000	NDIR Detector, Torch and Fusion
14-8318-150	3-port Series 66 Valve, 24 VDC
14-8907-350	Mass Flow Controller, THI, 1/8" Out
14-9593-090	Power Supply, 24 VDC
14-9639-090	Power Supply 25W
14-9640-090	Power Supply, 65 W, 12VDC
14-8781-290	DC Output Board
14-9180-190	PCB, Processor, Motor Fusion
14-9742-090	Display, LCD 2 x 20
14-9266-100	Carousel Drive Assy, Fusion/Torch
14-9624-000	Fusion/Torch Arm
14-9577-090	PCB, LED Banner Control Assembly
14-9298-028	Switch Power On/Off
14-7682-051	Solid State Relay, 240 VAC, 25A, 0 Volt

14-6273-000	Module, Power Entry, 250V 10A, Fil, fused
14-5505-034	Fuse, 2.5A, 250V, 5x20 mm, Slow Blow
14-9766-000	Keyboard / Mouse Combination
14-9729-091	Touchscreen Monitor Display (TOC Fusion)
14-9730-000	Fusion Internal PC

### Spare Parts - Tubing and Liquid / Gas Flow Path

14-9642-002	Tubing, Wire Reinforced (Clear), ea= 1ft
14-0440-002	Tubing, 1/16" OD, .04ID, NATL, FEP VAC BAKE
14-7630-002	Tubing, 1/8" OD, Red, Translucent (per inch)
14-7634-002	Tubing, 1/8" OD, Blue, Translucent (per inch)
14-7635-002	Tube, 1/8 OD, Green, Translucent
14-7633-002	Tubing, 1/8" OD, Yellow, Translucent (per inch)
14-7671-016	Ferrule, Super Flangeless, 1/16", PEEK
14-7695-016	Nut, PEEK, Natural for 1/16" OD Tubing
14-7038-016	Bushing, Exit, 1/8" Tefzel, Short
14-9778-016	Nut 1/8", 1/4 - 28, PEEK, Long
14-8659-016	Nut, Flangless, Short, 1/8, Peek
080-974	Union reducing 1/8 - 1/4 poly
14-9883-016	Union, Reducing, Poly Modified
14-2261-016	Union, 1/4 -1/16 Reducing, SS, W/ Teflon Fer
14-9596-053	Fusion, 3-Stage Needle, Septum Piercing
14-3228-001	Screw, Thumb, 6-32 X 3/8 LG
14-9620-094	Tray, Drain, Carousel
080-907	Tee 1/8 Polypropylene
080-974	Union reducing 1/8 - 1/4 poly
14-9751-067	Loop, MD Dryer, Teflon
090-099	Flow restrictor, 200CC +10 -5CC 1/8 Union Blu
14-9711-050	Pressure Regulator
14-9720-010	Nut Panel for Pressure Regulator
14-8661-001	Screw, 6-32 x 1-5 SHCS