

| Part Number | Description |
|---------------------------------------------------|---------------------------------------------------------------------------|
| Syringe, Valve, and Assemblies | |
| 14-9784-150 | Valve, 7 PORT (Phoenix and Apollo,) |
| 14-7514-001 | Syringe, 8 port valve plug |
| 14-7127-152 | Syringe, 25ml with plunger |
| 14-8587-30P | Syringe Pumper, Phoenix (NEW) |
| 14-8610-101 | Syringe plunger mounting screw, 8-32 shoulder screw, 5-port |
| 14-8438-179 | Syringe Stem |
| 14-3661-009 | #8x.09 Thick Washer |
| 110-252 | Knurled Thumb Screw 6-32 x 1/2 |
| 14-8640-080 | Bracket, Syringe Locking |
| 14-8661-001 | 6-32 x 1-5 Socket Head Cap Screw |
| 14-8505-086 | Encoder cable assembly, syringe drive pumper |
| 14-8455-047 | Timing Belt, syringe drive pumper |
| 14-8587-30pr | Syringe Pumper, Phoenix (Recertified) |
| Glassware and Assemblies in Chemistry Area | |
| 14-7015-024 | Reaction vessel, IC, for Phoenix 8000 |
| 14-7020-024 | Reaction vessel, UV, for Phoenix 8000 |
| 14-7021-100 | UV Reactor Cap Assembly (PEEK) |
| 14-7183-000 | UV-Lamp with amp connector and O-rings for Phoenix 8000 |
| 885-379 | Cap, sparger TEF |
| 14-7022-080 | Metal plate for UV reaction vessel cap |
| 14-7029-024 | Gas/liquid separator for Phoenix 8000 |
| 14-7033-043 | O-ring, sparge tube, in UV reaction vessel |
| 050-114 | Plug Con 6 CKT Amp 1480704-0 |
| 885-143 | Mist trap for Apollo 9000, DC-180 and Phoenix 8000 |
| 14-7014-024 | Scrubber assembly, chloride, for Apollo 9000 and Phoenix 8000 |
| 888-448 | Tube, scrubber U shape |
| 511-738 | Pyrex wool 1 gram |
| 511-876 | Tin, granular, 20 mesh - 1.5 oz. |
| 511-895 | Copper, granular, 20 mesh - 1.5 oz |
| 14-8203-000 | S03 Scrubber |
| 14-7680-024 | Fritted sparge tube, IC or UV, for Apollo 9000 or Phoenix 8000 |
| 100-446 | O-Ring, Lamp, UV Reaction Vessel |
| 110-236 | UV Reactor Screws (teflon cap) |
| 14-8809-000 | Tool for UV Lamp Lead Pins |
| 14-7193-000 | Shield, Tinted plastic cover with hinge (Smoked) |
| 110-206 | UV Reactor Screws (peek cap) |
| 14-3823-100 | Septum, with 1/8" hole, IC sparger |
| Tubing, Fittings, Septa, and Misc in | |
| 14-7637-002 | Tube, TFE 1/16" red translucent (price per inch) |
| 14-7681-002 | Tube, TFE 1/16" blue translucent (price per inch) |
| 14-0583-002 | Tubing, 1/8", teflon, FEP, (per foot) |
| 14-7630-002 | Tube, TFE 1/8" red translucent (price per inch) |
| 14-7635-002 | Tube, TFE 1/8" green translucent (price per inch) |
| 14-7634-002 | Tube, TFE 1/8" blue translucent (price per inch) |
| 14-7633-002 | Tube, TFE 1/8" yellow translucent (price per inch) |
| 154-424 | Tubing, Tygon 1/2" OD (523-359) (price per inch) |
| 517-799 | Septum, red/white for 1/16th in. tubing |
| 517-815 | Septa-connector, 1/32" hole, gray, ea (Formerly P/N 517-807, 10 per pack) |
| 517-813 | Septum, red for 1/8th in. tubing |
| 14-0063-043 | Septum, plug for mist trap |
| 14-8659-016 | Flangeless, short, 1/8", PEEK |
| 14-8820-016 | Flangeless, short, 1/16" |
| 14-7035-016 | Nut for 1/16" tubing |
| 14-7036-016 | Bushing, 1/16" for 8-port valve |
| 14-7013-094 | Drip tray for Phoenix 8000 |
| Detector | |
| 14-7674-2SL | IR, Non-linear, 220V/50Hz, 125mm cell |
| 010-962 | Cell liner, for Fuji and Apollo IR detectors, 125mm |

073-022 Cell window, for Fuji and Apollo IR detectors
 073-023 Cell window o-ring, for Fuji and Apollo IR detectors
 14-6860-000 Chopper motor/light source assembly for Fuji and Apollo IR detectors

Plumbing

14-0142-016 Connector, Clippard, 10-32-1/8"
 080-907 Tee 1/8 Polypropylene
 080-974 Union reducing 1/8 - 1/4 poly
 14-0356-016 Union, bulkhead, 1/8", filter assembly
 090-093 Flow Restrictor 20CC/Min Yellow
 090-097 Flow restrictor, 100mL/min, red dot
 090-099 Flow restrictor, 200 ml/min, blue dot
 14-7026-150 Valve, 2-way air, 24 VDC, for Phoenix 8000
 14-7025-150 Valve, 3-way air, 24 VDC, for Phoenix 8000
 14-7037-016 1/8" Flangless nut
 14-7038-016 Bushing, 1/8", Teflon
 090-825 Permeation dryer
 14-7195-000 Flowmeter assembly for Phoenix 8000
 14-4261-016 Cross union, 1/8", brass

Electronics

14-7009-090 Main CPU board for Apollo 9000 and Phoenix 8000
 14-7008-190 UV power supply for Phoenix 8000
 14-4383-028 Power entry module
 14-7007-090 Power supply (before S/N 00150001)
 14-7624-090 Power supply (after S/N 00150001)
 14-7717-100 Fan assembly
 14-7682-051 Solid state relay, 240 VAC, 25A, 0 Volt for Apollo, Phoenix, Fusion
 14-7735-000 UV transformer, 115V
 14-7735-100 UV transformer, 230V
 888-184 Solid state relay cover

Miscellaneous

14-7045-074 Phoenix 8000 operation manual
 14-8628-000 Filter Set / Maintenance Kit for Whatman 78-40 TOC Gas Generator (40-0014-000)
 14-9515-000 CO2 Filter for Regulator Assy, 14-9512-000 (Replace Annually)
 14-9512-000 Gas Regulator Assy - Includes Gas Regulator and CO2 Filter for carrier gas input into the Apollo.

DC-183 Boat Section

885-440 Phoenix & Apollo, PCBA 183 SSS, CTL.