

Item No. Description

Glassware

14-2333-024	Sparger, 25mL fritless, 1/2" (glassware only)
14-2334-024	Frit Sparger, 25mL, 1/2" neck, Tekmar
14-2336-024	Fritless Sparger, 5mL, 1/2" Tekmar
14-2337-024	Frit Sparger, 5mL, 1/2" Neck Tekmar
14-8924-100	25mL fritless Sparge Glassware Kit Atomx
14-8922-100	25mL frit Sparge Glassware Kit Atomx
15-0201-024	Vessel, 15mL Amber IS
15-0270-052	Syringe, 27 mL With Side Port

Heaters

14-9828-220	Trap Heater, Nichrome Wire, 230V, Atomx
14-9828-120	Trap Heater, Nichrome Wire, 115V Atomx
14-9895-120	Heater, Condensate Trap, 115V Atomx
14-9895-220	Heater, Condensate Trap, 230V, Atomx
14-8920-300	Assy, Heater, Sample 5/25mL 230VAC, Atomx
14-7806-000	Assy, RTD (24") Temp Sensor for Soil Valve, Mount
14-7806-100	Assy, RTD (28") Temp Sensor for Sample Heater
14-7806-200	Assy, RTD (34") Temp Sensor for Heated Water Reservoir
14-8105-320	Assy, Water Reservoir/Heater 110V
14-8785-100	Assy RTD, Oven, Atomx
14-8105-420	Assy, Water Reservoir/Heater 220V
14-8789-200	Assy, Mount HTR 115V Atomx
14-8790-E20	Assy, Transfer Line Heater 115V
14-8790-F20	Assy, Transfer Line Heater 230V
14-8920-500	Atomx, Sample Heater Kit (115V)
14-8920-600	Atomx, Sample Heater Kit (230V)
15-0100-100	Heater Assy, Sample Vial
15-0100-200	Heater Assy, Sample Vial, 55W, 230VAC
14-8789-400	Soil Valve Heater, Atomx 115V
14-8789-500	Soil Valve Heater, Atomx 230V
14-8788-300	Assy, Oven Heater, Atomx 115V

Fittings

12-0064-016	Bulkhead, 1/8" - 1/8" Tube stub, Brass
12-0070-016	Tee, 1/8", Brass
12-0075-016	Conn, Male, 1/8" Tube to 1/8" NPT
14-0051-016	Union, 1/16-1/16 SS
14-0241-016	Ferrule, 1/16", SS
14-0243-116	Nut, 1/16", Short Gold plated
14-1301-016	Ferrule, 1/2" Set, Teflon
14-1590-116	Nut, Male Plug, 1/16" Gold Plated
14-2261-016	Union, 1/4 - 1/16 Reducing, SS, W/ Teflon Fer
14-2535-016	Fitting, Tube 1/8" NPT 1/4", Brass
14-2792-016	Nut, Plug, 1/16", Brass
14-2950-116	Nut, Valco, 1" Long Gold Plated
14-3123-016	Ferrule, 1/16, Upchurch E/TFE (081-042)
14-3124-016	Nut, Short, Blk, 1/16 Upchurch
14-3354-016	Nut, 1/2" Knurled Swagelok
14-4050-016	Ferrule, 1/8 Fitt 1/8 Tube GVF
14-4261-016	Cross, 1/8" Union Brass Aps.
14-4602-116	Nut, 1/8", SS, Valco Gold Plated
14-4695-016	Tee, 1/16" Brass
14-5661-016	Pencil Filter, Filter Flow Module
14-7038-016	Bushing, Exit, 1/8" Tefzel, Short

14-7201-009	Washer, TEF
14-7514-001	Plug, Flat Point Nylon, 1/4-28 X 1/4" LG
14-7671-016	Ferrule, Super Flangeless, 1/16", PEEK
14-7695-016	Nut, PEEK, Natural for 1/16" OD Tubing
14-8470-016	Plug Nut, 1/4-28, Tefzel, Flat Bottom
14-8820-016	Nut, 1/16", Short, PEEK
14-8914-016	Nut, 1/8 Black Delrin Upchurch
14-9778-016	Nut 1/8", 1/4 - 28, PEEK, Long
14-9845-016	Trap Manifold Assy w/ Fittings
14-9924-016	Fitting, 1/4-28 to Barb
15-0088-016	Union, 1/4" Brass
15-0129-053	Needle, 3 Stage, Atomx
14-0158-016	Ferrule, 1/16", SS, Set
15-0171-016	Sample Mount Assy, w/o Fittings Atomx
15-0202-079	Bushing, Peek W/ O-Ring, IS Vessel
15-0273-000	Screen, 12 UM, Inline, IS Assembly

Electronics

14-5665-034	Fuse, Slow Blow, 5AMP, 250V, 5 X 20 MM
14-4383-028	Switch, Power, 10 AMP Filtered use Delta
14-5180-034	Fuse, 10A, 250V, 5x20mm, Slow Blow
14-8781-390	PCB, DC Control
14-9397-100	Optical Sensor Assy, 24"
14-9432-034	CAP Fuse, 2A
14-9433-034	CAP Fuse, 3.15A
14-9434-034	CAP Fuse, 5A
14-9577-090	PCB, LED Banner Control Assembly
14-9593-090	Power Supply, 24 VDC
14-9639-090	Power Supply 25W
14-9661-055	Encoder
14-9841-028	Switch, Limit, Miniature
15-0086-090	PCB, Temp Control
15-0197-090	PCB, Aux DC Board
15-0380-090	PCB, Motor Control

Valves

14-8317-100	Solid Valve Assy, Atomx
14-9719-000	Pumper, Kloehn, V6
14-9784-050	7-port valve for syringe pumper
14-9799-000	Mass Flow Controller Assy
15-0112-000	Manifold, Is Valve Atomx
15-0236-050	Valve, Is Atomx
15-0244-050	Regulator, Self relieving, 0-15psi
15-0271-000	Mounting Kit, IS Valve
15-0274-050	6-Port Valve Assy .030 Port W Conn
15-0275-250	6-Port Valve Actuator Assembly
15-0277-000	Clamp, Ring Valco Valve Heavy
15-0278-050	Check Valve (CV-3322)
15-0280-050	Valve, 24VDC, LQX 12

Miscellaneous

14-1197-039	Cord, Power, 115V 15 A, European Recept
14-8856-000	Assy, 10mL Antifoam 1520
14-5028-039	220V Power Cable
14-8834-500	Kit, Atomx Foam Sensor
14-8834-600	Kit, Atomx Eliminator Upgrade
15-0291-000	Chiller, Atomx
14-8920-500	Atomx, Sample Heater Kit (115V)

14-8920-600	Atomx, Sample Heater Kit (230V)
15-0000-017	Box and Foam, Atomx
15-0000-075	Atomx Teklink, CD
15-0000-KIT	Installation Kit Box
15-0279-000	Blank Water Reservoir, 10 Liters
15-0290-086	USB to RS232 Converter Cable
15-0293-000	Lubricant, Dupont Krytox #GPL206

Traps

12-0083-403	Trap, U Shape, Tenax 12" LG, Stamp "1"
12-0084-403	Trap, Tenax/SG, 12" LG, Stamp "2", U Shape
14-0124-403	Trap, Tenax/SG/Char 1/8" Stamp "3", U Shape
14-2366-403	Trap, Char-SC-TNX-OV1 12X1/8, "#5" U Shape
14-3347-403	Trap, OV-1. Tenax, #7 U Shape
14-3928-403	Strat-Trap, Carbo B/Carbo S, (#8)
14-5865-403	Strat-Trap, Vocarb 4000
14-5866-403	Strat-Trap, BTEX
14-9333-403	Strat-Trap, BTEX + MTBE
14-9908-403	Strat-Trap, #9 (U-shaped)
14-9909-403	Strat-Trap, Tenax/Silica Gel/Carbo Sieve III (#10)
14-5864-403	Strat-Trap, Vocarb 3000

Tubing

14-9927-002	Tubing, Siltek, 1/16" OD x 0.02" ID (per foot)
14-0583-002	Tubing, 1/8" OD, .062" ID, Natural FEP (per foot)
14-6934-002	PEEK Tubing, 1/16" OD, Natural (per foot)
14-7630-002	Tubing, 1/8" OD, Red, Translucent (per inch)
14-7633-002	Tubing, 1/8" OD, Yellow, Translucent (per inch)
14-7634-002	Tubing, 1/8" OD, Blue, Translucent (per inch)
14-8927-002	PEEK Tubing, 1/16" OD, Green (per foot)
14-9622-002	Tubing, Siltek, 1/16" OD x 0.040" ID (per foot)
15-0087-002	Tubing, Drain, Self Retracting (per foot)
15-0136-002	Tubing, 1/4" OD Blue Polyurethane (per foot)
14-8928-002	PEEK Tubing, 1/16" OD, Orange (per foot)
15-0272-002	Line Assy, I.S., w/ Minstac Connector

GC Interface Cables

14-2991-000	Interface, HP 5890 GC
14-2993-000	Interface, HP 5995/96/85/87/88/92 GC/MS with HP-1000/RTE GC/MS Software.
14-3147-000	Interface, general purpose/HNU 301/321/421 (Valve driver option necessary from HNU.)
14-3233-000	Interface, Perkin-Elmer 8000 Series/Autosystem, Clarus.
14-4188-086	Interface, HP 5890 w/5970 MSD and Unix or Pascal based software. Note: This cable requires the HP A111 (HP P/N 05990-60111) or A211 (HP P/N 05990-60211) accessory card and internal accessory cable (HP P/N 05987-60158) if using the Pascal software. If
14-4610-086	Interface, Shimadzu GC 14A/15A/17A GC 14A w/QP 1000 EX MSD and GC 14A w/QP 2000 MSD, also includes 1/16" to 2mm unions, needs Shimadzu P/N 221-34280-91
14-4610-186	Interface, Shimadzu GC 14A/15A/17A GC 14A w/QP 1000 EX MSD and GC 14A w/QP 2000 MSD, needs Shimadzu P/N 221-34280-91
14-4652-086	Interface, HP 5890 w/5970 MSD and Unix-B or MS-DOS software, HP 5890/5971/5972 MSD and Unix-B or MS-DOS software and HP 5890/5989 MS Engine
14-4938-086	Interface, Carlo Erba Mega and Vega Series/Fisons 8000

14-5397-086	Interface, two Tekmars on one Perkin-Elmer 8000 Series/Autosystem
14-6689-086	Interface, HP 6850, 6890, 7890 GC
14-7396-086	Interface, Varian 3800/3900
14-7397-086	Trace 2000/Trace GC Ultra Interface cable w/ 1/16" to 2mm unions
14-7397-186	Trace 2000/Trace GC Ultra Interface cable
14-8417-000	Duet™ High Throughput Interface. Interfaces two Teledyne Tekmar VOC systems with any GC or GC/MS system.
14-3233-200	PE Clarus 500 Interface Cable
14-9754-086	Shimadzu 2010-2014 Interface Cable

Cables

14-1197-039	Cord, Power, 115V 15 A, European Recept
14-5107-086	RS232 Cable Assembly (included with unit)
14-8002-186	Cable, Assembly, DC Power Harness
14-8797-386	Cable Assy, CPU to SPI Interface
14-8798-386	Cable Assy, GC/IO Interface
14-8876-386	Cable Assy, Mass Flow Controller Comm/Power

Mechanisms

14-9401-000	Stepper motor Assy
14-9571-022	Latch, Front Cover
-547817	Elevator Assy, Atomx
15-0076-280	Tray, Carousel, 80 Position Atomx
15-0110-000	Belt, Mixer Drive, Atomx
15-0180-300	Vial Retraction Assy Atomx
15-0185-079	Bearing Guide, Carousel Tray Atomx
15-0216-000	Assy, Elevator Sample Motor, Atomx
15-0233-031	Knob, Carousel Tray Atomx
15-0255-000	Carousel, Drive Assy, Atomx
15-0259-000	Guard, Assy Mount Plastic, Atomx