

Item No.	Description
Caps	
14-4436-000	Crimp Cap. For headspace septa and 20mm top headspace vials, 125 per pack.
Septa	
14-7418-080	Foil backing, septa (pack of 125), required for high temperature vials.
14-4385-043	Butyl rubber septa without teflon face. For 20mm top headspace vials, 125 per pack. Recommended for temperatures up to 100 degrees C. Color coded grey. Silicon rubber/teflon face septa. For 20mm top headspace vials, 125 per pack. Included with instrument shipment. (Recommended for temperatures over 100 degrees C. Color coded tan.)
14-5818-043	
Vials	
14-4440-324	Vials, 22ml headspace, w/20mm top, for high temperature use (over 225C), 125 per pack. Included with instrument shipment.
14-4440-024	Vials, 22ml headspace, w/20mm top, 125 per pack.
14-4439-024	Vials, 12ml headspace, w/20mm top, 125 per pack. ** Requires insert for use**
14-4438-024	Vials, 9ml headspace, w/20mm top, 125 per pack. ** Requires insert for use**
Accessories	
14-9133-079	Vial insert for 9ml and 12ml vials. (1 insert required per carousel position)
14-9133-179	Vial insert for 11ml vials. 10pk (1 insert required per carousel position)
Electronics	
14-1197-039	Cord, Power, 115V 15 A, European Recept
14-5107-086	RS232 Cable Assembly (included with unit)
14-8313-190	Power Supply, 24 VDC, 300W-Solatek Only
14-8312-090	Power Supply, triple out, +5/+15/-15 VDC
14-7835-090	Temperature Control Board
14-9185-090	PCBA, AC Output Board
14-9180-090	PCBA, Processor, Motor
14-8002-186	Cable, Assembly, DC Power Harness
14-8797-086	Cable, Assembly, CPU to SPI Interface.
14-8798-186	Cable, Assembly, GCIO Interface.
14-8876-186	Cable, Assembly, MFC Comm.
14-9396-086	Cable, Assembly, Valco/MFC Power.
14-9181-086	Cable, Assembly, CPU to Display Cable.
14-8003-186	Cable, Assembly, 24 VDC.
14-8004-386	Cable, Assembly, 6 lead Telephone line.
14-9403-086	Cable, Assembly, Main AC Harness.
14-8020-086	Temp to Solid State relay board cable
14-9407-086	Cable, Assembly, Valco Control.
14-8876-086	Mass Flow Controller Cable
14-7105-086	Cable, host, RS 232
14-8904-586	Assembly, valve LED power cable
14-9437-000	Cooling fan, front and rear
14-9477-000	RTD Jumper
14-9404-100	Platen RTD Assy
14-9404-000	Oven RTD Assy, High Temp Leads
14-8191-000	Assy, lighted display
14-4383-028	Switch, Power, 10 AMP Filtered use Delta
14-5180-034	Fuse, 10A, 250V, 5x20mm, Slow Blow
14-5665-034	Fuse, Slow Blow, 5AMP, 250V, 5 X 20 MM
14-9139-000	Dual Optical Sensor Assy

14-8781-290 DC Output Board

Fittings

14-0243-016	Nut, 1/16", short, SS, Valco
14-0241-016	Ferrule, 1/16", SS
14-8916-016	Nut, Black ,Upchurch, 1/16 Long
14-3123-016	Ferrule, 1/16, Upchurch ETFE (081-042)
14-9410-016	Quick Disconnect, 1/16".
14-0264-016	Union, bulkhead, 1/16", SS, Upchurch
14-1590-116	Nut, Male Plug, 1/16" Gold Plated
12-0070-016	Tee, 1/8", Brass
14-0051-016	Union, 1/16-1/16 SS
14-0158-016	Ferrule, 1/16", SS, Set
14-3295-016	Nut, 1/16", male, SS, Swagelok
14-2792-016	Nut, Plug, 1/16", Brass
14-3296-016	Nut, 1/16", SS, Swagelok, plug
14-8470-016	Plug Nut, 1/4-28, Tefzel, Flat Bottom
14-0243-116	Nut, 1/16", Short Gold plated
14-7671-016	Ferrule, Super Flangeless, 1/16", PEEK
14-7695-016	Nut, PEEK, Natural for 1/16" OD Tubing

Instrument Interface Options

14-6689-086	Interface, HP 6850, 6890, 7890 GC
14-2991-000	Interface, HP 5890 GC Interface, HP 5890 w/5970 MSD and Unix or Pascal based software. Note: This cable requires the HP A111 (HP P/N 05990-60111) or A211 (HP P/N 05990-60211) accessory card and internal accessory cable (HP P/N 05987- 60158) if using the Pascal software. If using Unix software the internal accessory cable (HP P/N 05987-60158) only is necessary.
14-4188-086	Interface, HP 5890 w/5970 MSD and Unix-B or MS-DOS software, HP 5890/5971/5972 MSD and Unix-B or MS-DOS software and HP 5890/5989 MS Engine
14-4652-086	Interface, HP 5995/96/85/87/88/92 GC/MS with HP-1000/RTE GC/MS Software.
14-2993-000	Interface, Perkin-Elmer 8000 Series/Autosystem, Clarus.
14-3233-000	Interface, Varian 3800/3900
14-7396-086	Interface, Shimadzu GC 14A/15A/17A GC 14A w/QP 1000 EX MSD and GC 14A w/QP 2000 MSD, also includes 1/16" to 2mm unions, needs Shimadzu P/N 221-34280-91
14-4610-086	Interface, Shimadzu GC 14A/15A/17A GC 14A w/QP 1000 EX MSD and GC 14A w/QP 2000 MSD, needs Shimadzu P/N 221-34280-91
14-4610-186	Interface, Carlo Erba Mega and Vega Series/Fisons 8000
14-4938-086	Interface, general purpose/HNU 301/321/421 (Valve driver option necessary from HNU.)
14-3147-000	Trace 2000/Trace GC Ultra Interface cable w/ 1/16" to 2mm unions
14-7397-086	Trace 2000/Trace GC Ultra Interface cable
14-7397-186	Trace 2000/Trace GC Ultra Interface cable
14-3233-200	PE Clarus 500 Interface Cable
14-9754-086	Shimadzu 2010-2014 Interface Cable
14-7837-690	CPU Board (for use with Stratum)

Loops

14-9172-067	HT3 1ml loop, Silco coated.
14-9417-067	HT3 2 ml loop, Silco coated.
14-8641-067	Loop, 5ml, Silcosteel
14-9414-067	HT3 100ul loop, Silco coated.
14-9416-067	HT3 500ul loop, Silco coated.
14-9415-167	HT3 300ul loop, Silco coated
14-9415-067	HT3 250ul loop, Silco coated.

Sample Preparation Accessories

14-4863-027	Hand crimper for 20mm top headspace vials
14-5020-027	Crimp station, air powered
14-4864-027	Tool, decapper 20mm neck
14-4865-027	Crimp station, 20mm top headspace vials. Includes 20mm jaw set
14-4386-027	Jaw set, 20mm, crimp station
14-4866-027	Rack holder, vial, for 10 headspace vials
14-9001-077	Electronic Crimper for 12x32 mm vials.(115V)
14-9001-177	Electronic Crimper for 12x32mm vials. (230V)

TekLink HT3

14-9300-075	HT3 TekLink Software.
14-5745-090	Serial port expansion board for PC (required if your computer has no available serial ports)

Tubing

14-8927-002	PEEK Tubing, 1/16" OD, Green (per foot)
14-7634-002	Tubing, 1/8" OD, Blue, Translucent (per inch)
14-8938-002	Tubing, Silcosteel, 1/16"OD x .04" ID (no siloxane)
14-8928-002	PEEK Tubing, 1/16" OD, Orange (per foot)
14-6934-002	PEEK Tubing, 1/16" OD, Natural (per foot)
14-9384-002	PEEK Tubing, restrictive, 0.010 ID x 1/16" OD, Blue (per foot)

Valves and Pneumatics

14-8318-150	3-port Series 66 Valve, 24 VDC
14-9145-050	8-Port Valve Assembly
14-8907-000	Mass Flow Controller Assy, HT3
14-9145-150	8-port valve only
14-9145-250	8-port valve rotor
14-9145-350	8-port valve and actuator
14-9145-450	8-port valve control box
14-9123-000	Needle stripper assy

Computer Interface Cables

14-5107-086	RS232 Cable Assembly (included with unit)
14-9182-086	Cable Assy, RJ45 Ethernet

Heaters

14-8791-100	Trap Heater assembly, 110V
14-8791-200	Trap Heater assembly, 230V
14-8790-920	Transfer line assy, 115V
14-8790-820	Transfer line assy; 230V
14-9183-000	Oven heater cartridge assy, 115V
14-9183-100	Oven heater cartridge assy, 230V
14-9085-120	TRL heater assy, 115V
14-9085-220	TRL heater assy, 230V
14-9240-100	Platen heater assembly, 115 VAC
14-9240-200	Platen heater assembly, 230V

Motor Assemblies

14-9131-000	Vial loader motor assy
14-9155-000	Sample/mixer motor assy
14-9266-000	Carrousel drive assy
14-9258-000	Platen worm assembly
14-9401-000	Stepper motor Assy